

RYS/NLC-L SPO2 EAR SENSOR WITHOUT ADAPTER

SunTech Medical®

Réf. DSU250906



RYS/NLC-L SPO2 EAR SENSOR WITHOUT ADAPTER

Réf. **DSU250906**

FEATURES			
Medical device	No		
Marquage		Barcode	3661040217699
Guarantee (month):	24	Customs Code	90181990
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	2,0
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	16,0
Weight (kg) :		Weight packaging (kg) :	0,010
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:47 - This article is available on our website : www.gsh-med.fr with the following reference DSU250906

All characteristics are given for information only and may be changed