

## PROFESSIONAL PERSONAL FLOOR SCALE 300KG - CLASS III


Réf. EK MPC02-CIII

### KERN

GARANTIE  
**2**  
ANS  
YEARS  
GUARANTEE

CE  
0297

DISPOSITIF MEDICAL  
MEDICAL DEVICE  
CLASSE  
**I<sub>m</sub>**



# PROFESSIONAL PERSONAL FLOOR SCALE 300KG - CLASS III



Réf. EKMPC02-CIII

**Version 300 kg and larger weighing plate! Specially suitable for weighing overweight patients - Class III**

The bestseller for doctor's surgeries, chemists and hospitals

- Robust design for daily use in a professional environment
- Ergonomically optimised design
- For easy and hygienic cleaning
- Secure and non-slip positioning with height-adjustable rubber feet
- Level indicator to level the balance precisely
- KERN MPC 300K-1LM: Sturdy handle for easy transportation of the scale
- Flexible positioning of display device: e.g. free-standing or mounted to the wall. Wall mount for display device, standard
- Large platform with non-slip and wear-resistant surface made of plastic, to provide a secure footing for patients
- Hold function: When patients do not stand completely still, a stable weight is calculated using an average weight and this is then "frozen". This means that you have sufficient time to attend to the patient first and then take the weight reading in peace
- BMI function to determine under weight/normal weight/over weight
- Mother and child function: the mother's weight will be taken and then the scale will be reset to zero at the push of a button. Then the mother and baby can be weighed together. The weight of the child will now be shown on the display. This is particularly good for infants, as this way they can be cradled, safe and sound, in the arms of their mother while being weighed
- Battery- or mains-powered, Optional rechargeable battery pack
- Reference for same scale without Verification class III: EKMPC02

## FEATURES

Medical device	Yes	Classe Im (fonction de mesurage)	
Marquage	CE0297	Barcode	3661040217712
Guarantee (month):	24	Customs Code	
Sterile	No		

## INFORMATIONS

Main unit	Packaged unit
Width (cm) :	Width packaging (cm):
Length (cm) :	Length packaging (cm) :
Height (cm) :	Height packaging (cm) :
Weight (kg) :	Weight packaging (kg) :

## SPECIFIC STORAGE CONDITION

Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference EKMPC02-CIII

All characteristics are given for information only and may be changed

## PROFESSIONAL PERSONAL FLOOR SCALE 300KG - CLASS III



Réf. EKMPC02-CIII

### COMPOSITION

Code	Descriptive	Quantity
EKCERTI06	EU-Verification (CE) electr.balances 50 kg - 350 kg	1
EKMPC02	Personal floor scale MPC 300kg, 100gr.	1