

BATTERY FOR ES 60 / PRO M / PRO H SPENGLER



Réf. 550616



BATTERY FOR ES 60 / PRO M / PRO H SPENGLER



Réf. 550616

FEATURES			
Medical device	No		
Marquage		Barcode	3700446008901
Guarantee (month):	0	Customs Code	85075000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	5,7	Width packaging (cm):	
Length (cm) :	5,1	Length packaging (cm) :	
Height (cm) :	2,8	Height packaging (cm) :	
Weight (kg) :	0,119	Weight packaging (kg) :	0,120
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:52 - This article is available on our website : www.gsh-med.fr with the following reference 550616

All characteristics are given for information only and may be changed