

**POWER SUPPLY (CH,EURO,UK) 12V, 500MA,
INPUT 100V-240V (FROM 2017)**

KERN

Réf. EKTRAN01

POWER SUPPLY (CH,EURO,UK) 12V, 500MA, INPUT 100V-240V (FROM 2017)



Réf. **EKTRAN01**

FEATURES			
Medical device	No		
Marquage		Barcode	3661041009925
Guarantee (month):	0	Customs Code	85044090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,0	Width packaging (cm):	
Length (cm) :	10,0	Length packaging (cm) :	
Height (cm) :	6,0	Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:14 - This article is available on our website : www.gsh-med.fr with the following reference EKTRAN01

All characteristics are given for information only and may be changed