

## MASTER PALM 2 - HANDHELD OXIMETER W/ INTERCHANGEABLE SENSORS + COVERS



Réf. 223000P



## MASTER PALM 2 - HANDHELD OXIMETER W/ INTERCHANGEABLE SENSORS + COVERS



Réf. **223000P**

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0123	Barcode	
Guarantee (month):	24	Customs Code	9018199000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,718	Weight packaging (kg) :	0,718
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min	15	Humidity max	85



07/04/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 223000P

All characteristics are given for information only and may be changed

## MASTER PALM 2 - HANDHELD OXIMETER W/ INTERCHANGEABLE SENSORS + COVERS



Réf. 223000P

### COMPOSITION

Code	Descriptive	Quantity
223000	MASTER PALM 2 - Handheld Oximeter With Interchangeable Sensors	1
223004	Silicon Protect Cover for MASTER PALM 2/3 (2 pcs)	1