


SOFTWARE PC-ECG FOR ECG CARDIOMATE (WITH PHYSICAL KEY)



Réf. 260038



SOFTWARE PC-ECG FOR ECG CARDIOMATE (WITH PHYSICAL KEY)



Réf. **260038**

FEATURES			
Medical device	No		
Marquage		Barcode	3700446042073
Guarantee (month):	12	Customs Code	9018110000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



08/07/2022 12:44 - This article is available on our website : www.gsh-med.fr with the following reference 260038

All characteristics are given for information only and may be changed