

Réf. DTBAND01



CE

DMDIV

Réf. DTBAND01

FEATURES			
Medical device	Yes	DMDIV (diagnostic In Vitro)	
Marquage	CE	Barcode	3760309510132
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	6,0
Length (cm) :		Length packaging (cm) :	16,0
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :	0,100	Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)


 31/03/2022 03:39 - This article is available on our website : www.gsh-med.fr with the following reference DTBAND01

All characteristics are given for information only and may be changed