

Réf. DSU230701



Réf. DSU230701

| FEATURES | | | |
|----------------------------|-----------|-------------------------|----------------|
| Medical device | Yes | Classe IIa | |
| Marquage | CE0413 | Barcode | 10840935105718 |
| Guarantee (month): | 24 | Customs Code | 9018.90.5040 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 17,4 | Width packaging (cm): | |
| Length (cm) : | 24,0 | Length packaging (cm) : | |
| Height (cm) : | 11,5 | Height packaging (cm) : | |
| Weight (kg) : | 1,680 | Weight packaging (kg) : | |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 03:41 - This article is available on our website : www.gsh-med.fr with the following reference DSU230701

All characteristics are given for information only and may be changed