

SUNTECH CT40 PNI + SPO2 NELLCOR COMPATIBLE + TEMPERATURE COVIDIEN FIL

SunTech Medical®

Réf. DSU551443N



GARANTIE
3
ANS
YEARS
GARANTEE

CE
0413

DISPOSITIF MEDICAL CLASSE II
MEDICAL DEVICE

SUNTECH CT40 PNI + SPO2 NELLCOR COMPATIBLE + TEMPERATURE COVIDIEN FIL

Réf. DSU551443N

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0413	Barcode	
Guarantee (month):	36	Customs Code	9018199560
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:44 - This article is available on our website : www.gsh-med.fr with the following reference DSU551443N

SUNTECH CT40 PNI + SPO2 NELLCOR COMPATIBLE + TEMPERATURE COVIDIEN FIL

Réf. DSU551443N

COMPOSITION

Code	Descriptive	Quantity
DSU551400		1
DSU551514	DELUXE MOBILE STAND CT40	1
DSU551401		1
DSU551403		1