

## PE/PLASTIC SCOOP STRETCHER, YELLOW

Réf. SBRAN13\_J



## PE/PLASTIC SCOOP STRETCHER, YELLOW



Réf. SBRAN13\_J

- «Spoon» function performed by the two longitudinal blades, to position the stretcher under the patient without moving him. The scoop stretcher can be used to lift, move, or transport patients.
- Length adjustable on 4 levels according to the patient's size, from 162 to 234 cm
- High resistance aluminium structure
- Dimensions (L x W x H): 162-234 x 44 x 7 cm
- Weight: 7 kg
- Maximum load carried: 220 kg

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040094670
Guarantee (month):	12	Customs Code	94029000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	45,0	Width packaging (cm):	47,0
Length (cm) :	162,0	Length packaging (cm) :	165,0
Height (cm) :	7,0	Height packaging (cm) :	8,0
Weight (kg) :	8,000	Weight packaging (kg) :	10,800
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:21 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SBRAN13\_J