

ALUMINIUM ALLOY STAIR STRETCHER -W 38 CM - 2 WHEELS, ORANGE

Réf. SBRAN14_O



ALUMINIUM ALLOY STAIR STRETCHER -W 38 CM - 2 WHEELS, ORANGE



Réf. **SBRAN14_O**

- High strength aluminium alloy structure
- Ideal for transferring a patient, climbing up and down the stairs
- Two fixed rear wheels allow movement on most surfaces with minimal effort
- The safety of patients during the transport process is ensured by 3 restraining straps with quick release buckle
- At the front of the stretcher, two retractable arms
- Dimensions (L x W x H): 62 X 51 x 91 cm, unfolded; 90 x 17 x 59 cm, folded
- Weight: 8 kg
- Maximum load carried: 159 kg
- Soft seat
- Seat height: 49 cm
- Single colour: orange

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe I | |
| Marquage | CE | Barcode | 3661040094687 |
| Guarantee (month): | 12 | Customs Code | 94029000 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | 52,0 | Width packaging (cm): | 56,0 |
| Length (cm) : | 84,0 | Length packaging (cm) : | 98,0 |
| Height (cm) : | 92,0 | Height packaging (cm) : | 16,0 |
| Weight (kg) : | 9,000 | Weight packaging (kg) : | 10,000 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 23:21 - This article is available on our website : www.gsh-med.fr with the following reference SBRAN14_O