

ARRAL TROUSERS, WHITE, SIZE 00



Réf. WARR07_00



ARRAL TROUSERS, WHITE, SIZE 00



Réf. WARR07_00

FEATURES			
Medical device	No		
Marquage	Barcode	3661040033488	
Guarantee (month):	12	Customs Code	62034311095
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	28,0	Width packaging (cm):	28,0
Length (cm) :	32,0	Length packaging (cm) :	32,0
Height (cm) :	1,0	Height packaging (cm) :	1,0
Weight (kg) :	0,330	Weight packaging (kg) :	0,330
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:21 - This article is available on our website : www.gsh-med.fr with the following reference WARR07_00