

ARRAL TROUSERS, WHITE, SIZE 5



Réf. WARR07_95



ARRAL TROUSERS, WHITE, SIZE 5



Réf. WARR07_95

FEATURES			
Medical device	No		
Marquage		Barcode	3661040033549
Guarantee (month):	12	Customs Code	62034311095
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	37,0	Width packaging (cm):	37,0
Length (cm) :	42,5	Length packaging (cm) :	42,5
Height (cm) :	2,0	Height packaging (cm) :	2,0
Weight (kg) :	0,330	Weight packaging (kg) :	0,330
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 23:21 - This article is available on our website : www.gsh-med.fr with the following reference WARR07_95

All characteristics are given for information only and may be changed

2 / 2

HOLTEx

30 rue Jean de Guiramand
13290 Aix-en-Provence - France

www.gsh-med.fr

Contact :

+33(0)4 42 90 31 31 - info@gsh-med.fr