

PRE-CONNECTED ELECTRODE

SaverOne®

Réf. SJACCE01



CE

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
I



PRE-CONNECTED ELECTRODE



Réf. SJACCE01

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040055145
Guarantee (month):	0	Customs Code	90189085
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	17,0
Length (cm) :	21,0	Length packaging (cm) :	20,0
Height (cm) :	1,0	Height packaging (cm) :	3,0
Weight (kg) :	0,128	Weight packaging (kg) :	0,128
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:12 - This article is available on our website : www.gsh-med.fr with the following reference SJACCE01

All characteristics are given for information only and may be changed