

## PEDIATRIC ELECTRODE

SaverOne®

Réf. SJACCE03



CE

DISPOSITIF MEDICAL  
MEDICAL DEVICE

CLASSE  
I

## PEDIATRIC ELECTRODE



Réf. **SJACCE03**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040055169
Guarantee (month):	0	Customs Code	90189085
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	17,5	Width packaging (cm):	17,0
Length (cm) :	25,5	Length packaging (cm) :	23,0
Height (cm) :	5,0	Height packaging (cm) :	3,0
Weight (kg) :	0,120	Weight packaging (kg) :	0,120
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:13 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SJACCE03

All characteristics are given for information only and may be changed