


PUNCH BIOPSY Ø 3 MM BOX OF 10

Réf. CPPUNC05



PUNCH BIOPSY Ø 3 MM BOX OF 10



Réf. **CPPUNC05**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE2460	Barcode	3661040186191
Guarantee (month):	12	Customs Code	90183290
Sterile	Yes		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :	8,2	Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	3,5
Weight (kg) :		Weight packaging (kg) :	0,064
SPECIFIC STORAGE CONDITION			
Temperature min	10	Temperature max	40
Humidity min	35	Humidity max	65



07/04/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference CPPUNC05