


PUNCH BIOPSY Ø 2 MM BOX OF 10

Réf. CPPUNC03



PUNCH BIOPSY Ø 2 MM BOX OF 10



Réf. **CPPUNC03**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE2460	Barcode	3661040186177
Guarantee (month):	12	Customs Code	90183290
Sterile	Yes		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	8,2	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,064
SPECIFIC STORAGE CONDITION			
Temperature min	10	Temperature max	40
Humidity min	35	Humidity max	65



07/04/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference CPPUNC03