

SILICONE TUBE ROLL OF Ø6X10 MM - 1 METER

Réf. BCASPI201



CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEED

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE **II** CLASSE

SILICONE TUBE ROLL OF Ø6X10 MM - 1 METER



Réf. BCASPI201

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040166353
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,0	Width packaging (cm):	
Length (cm) :	100,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:12 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI201