

SILICONE TUBE ROLL OF Ø8X14 MM - 1 METER

Réf. BCASPI202

CA-MI



CE
0123



SILICONE TUBE ROLL OF Ø8X14 MM - 1 METER



Réf. BCASPI202

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	3661040169132
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:12 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI202

2/2