

## S/S BOX 30 X 12 X 6 CM

Réf. AI30126



## S/S BOX 30 X 12 X 6 CM



Réf. AI30126

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040019598
Guarantee (month):	12	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,0	Width packaging (cm):	12,0
Length (cm) :	30,0	Length packaging (cm) :	30,0
Height (cm) :	6,0	Height packaging (cm) :	6,0
Weight (kg) :	0,758	Weight packaging (kg) :	0,810
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 14:48 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference AI30126

All characteristics are given for information only and may be changed