

NEW ASKIR 20 ELECTRICAL SUCTION UNIT, 1L AUTOCLAVABLE POLYCARBONATE J

Réf. BCASPI06

CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL CLASSE II
MEDICAL DEVICE



NEW ASKIR 20 ELECTRICAL SUCTION UNIT, 1L AUTOCLAVABLE POLYCARBONATE J

CA-MI

Réf. BCASPI06

FEATURES

Medical device Yes Classe IIa

Marquage CE0123

Barcode 8033622540137

Guarantee (month): 24

Customs Code 90189084

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) : 35,0

Width packaging (cm):

Length (cm) : 21,0

Length packaging (cm) :

Height (cm) : 18,0

Height packaging (cm) :

Weight (kg) : 2,500

Weight packaging (kg) : 3,550

SPECIFIC STORAGE CONDITION

Temperature min -20

Temperature max 45

Humidity min 1

Humidity max 85



30/03/2022 22:43 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI06

All characteristics are given for information only and may be changed

2 / 2

HOLTEX

30 rue Jean de Guiramand
13290 Aix-en-Provence - France

www.gsh-med.fr

Contact :

+33(0)4 42 90 31 31 - info@gsh-med.fr