

**CONTAINER FLOVAC 2L (REUSABLE) FOR  
ASPIRET, ASKIR 20, 30 AND C30BR**



Réf. BCASPI104

# CONTAINER FLOVAC 2L (REUSABLE) FOR ASPIRET, ASKIR 20, 30 AND C30BR



Réf. **BCASPI104**

FEATURES			
Medical device	No		
Marquage		Barcode	0000000000611
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	34,0
Height (cm) :		Height packaging (cm) :	14,0
Weight (kg) :		Weight packaging (kg) :	0,366
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:10 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI104