

SET OF SILICONE TUBES Ø 6X10 FLOVAC + CONNECTOR ASPIRET, ASKIR 20/30

CA-MI

Réf. BCASPI111

GARANTIE
1
AN
TVA
GARANTIE



CE
0123

SET OF SILICONE TUBES Ø 6X10 FLOVAC + CONNECTOR ASPIRET, ASKIR 20/30



Réf. BCASPI111

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	3661040160818
Guarantee (month):	12	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,5	Width packaging (cm):	
Length (cm) :	11,0	Length packaging (cm) :	
Height (cm) :	3,5	Height packaging (cm) :	
Weight (kg) :	0,050	Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:09 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI111

All characteristics are given for information only and may be changed