

## SILICONE TUBE ROLLER Ø6X10MMX50M NEW ASPIRET / ASKIR / EMIVAC

Réf. BCASPI144

CA-MI



CE  
0123



# SILICONE TUBE ROLLER Ø6X10MMX50M NEW ASPIRET / ASKIR / EMIVAC

Réf. BCASPI144

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	0000000000895
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	41,0
Length (cm) :		Length packaging (cm) :	41,0
Height (cm) :		Height packaging (cm) :	13,0
Weight (kg) :		Weight packaging (kg) :	3,500
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI144