

ASKIR 118 WM (WALL MOUNTED) SUCTION UNIT, 1L JAR IN MAKROLON CA-MI

Réf. BCASPI306

CA-MI

GARANTIE
3
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II
a



ASKIR 118 WM (WALL MOUNTED) SUCTION UNIT, 1L JAR IN MAKROLON CA-MI



Réf. **BCASPI306**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	8054619100646
Guarantee (month):	36	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	35,0	Width packaging (cm):	
Length (cm) :	15,0	Length packaging (cm) :	
Height (cm) :	19,0	Height packaging (cm) :	
Weight (kg) :	2,700	Weight packaging (kg) :	5,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:31 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI306