

DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR

Réf. BCASPI101

CA-MI

GARANTIE
1
AN
TIAS
GARANTIE



DISPOSABLE FLOVAC LINER 2L FOR ASPIRET, ASKIR 20, 30, C30BR



Réf. BCASPI101

FEATURES			
Medical device	No		
Marquage	Barcode	0000000000604	
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:10 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI101