

12V CAR ADAPTER FOR ASKIR BR

CA-MI

Réf. BCASPI98

12V CAR ADAPTER FOR ASKIR BR

CA-MI

Réf. BCASPI98

FEATURES			
Medical device	No		
Marquage		Barcode	3661040169446
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	8,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	0,076
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	5	Humidity max	85



30/03/2022 23:37 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI98