

## SET OF SILICONE TUBES Ø8X14 JAR + FILTER FOR ASKIR C30BR

Réf. BCASPI116

CA-MI



CE  
0123



# SET OF SILICONE TUBES Ø8X14 JAR + FILTER FOR ASKIR C30BR

Réf. **BCASPI116**

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	3661040166346
Guarantee (month):	12	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,242
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



23/09/2022 21:07 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI116