

## AUTOCLAVABLE 4000 / 5000ML JAR HOLDER NEW HOSPIVAC

Réf. BCASPI132

CA-MI

GARANTIE  
1  
AN  
TVA  
GARANTIE



## AUTOCLAVABLE 4000 / 5000ML JAR HOLDER NEW HOSPIVAC

CA-MI

Réf. BCASPI132

FEATURES			
Medical device	No		
Marquage		Barcode	0000000000802
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	20,0
Length (cm) :		Length packaging (cm) :	40,0
Height (cm) :		Height packaging (cm) :	6,0
Weight (kg) :		Weight packaging (kg) :	0,388
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI132

All characteristics are given for information only and may be changed

2/2