

## ADAPTER GASKET FOR 2000ML JAR NEW HOSPIVAC

Réf. BCASPI133



## ADAPTER GASKET FOR 2000ML JAR NEW HOSPIVAC

CA-MI

Réf. BCASPI133

FEATURES			
Medical device	No		
Marquage		Barcode	0000000000819
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI133

2/2