

ROLL OF SILICONE TUBE Ø6X10MM LENGTH 10M FOR ASPIRET, ASKIR

Réf. BCASPI143

CA-MI



CE
0123



ROLL OF SILICONE TUBE Ø6X10MM LENGTH 10M FOR ASPIRET, ASKIR



Réf. **BCASPI143**

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	0000000000888
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	21,0
Length (cm) :		Length packaging (cm) :	23,0
Height (cm) :		Height packaging (cm) :	5,0
Weight (kg) :		Weight packaging (kg) :	0,694
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	10	Humidity max	85



30/03/2022 23:36 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI143