

SET OF SILICONE TUBES FOR AUTOCLAVABLE JARS FOR NEW ASPIRET, ASKIR

Réf. BCASPI136

CA-MI



CE
0123



SET OF SILICONE TUBES FOR AUTOCLAVABLE JARS FOR NEW ASPIRET, ASKIR

CA-MI

Réf. BCASPI136

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	3661040163604
Guarantee (month):	12	Customs Code	90330090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	4,340
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:36 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI136

All characteristics are given for information only and may be changed

2/2