

NEW ASPIRET ELECTRICAL SUCTION UNIT 1L AUTOCLAVABLE POLICARBONATE JAR

Réf. BCASPI01



CA-MI

GARANTIE
2
ANS
YEARS
GUARANTEE

CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE
CLASSE
II
B

NEW ASPIRET ELECTRICAL SUCTION UNIT 1L AUTOCLAVABLE POLICARBONATE JAR



Réf. **BCASPI01**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	8033622540229
Guarantee (month):	24	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	23,0	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	16,0	Height packaging (cm) :	
Weight (kg) :	2,200	Weight packaging (kg) :	2,800
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 22:43 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI01