

## AMBULANCE WALL MOUNT FOR ASKIR 118VM

Réf. BCASPI169

CA-MI



CE  
0123



## AMBULANCE WALL MOUNT FOR ASKIR 118VM



Réf. BCASPI169

FEATURES			
Medical device	No		
Marquage	CE0123	Barcode	3661040169118
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	20,5
Length (cm) :		Length packaging (cm) :	25,0
Height (cm) :		Height packaging (cm) :	16,0
Weight (kg) :		Weight packaging (kg) :	0,802
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	1	Humidity max	85



30/03/2022 23:36 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference BCASPI169