

VACUUM GAUGE FOR ASPIRET

Réf. BCASPI149



VACUUM GAUGE FOR ASPIRET



Réf. **BCASPI149**

FEATURES			
Medical device	No		
Marquage		Barcode	0000000000932
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	5,5
Length (cm) :		Length packaging (cm) :	5,0
Height (cm) :		Height packaging (cm) :	5,5
Weight (kg) :		Weight packaging (kg) :	0,044
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	45
Humidity min	0	Humidity max	85



30/03/2022 23:36 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI149