

SILICONE TUBES 8X14MMX150CM/8X14MMX17+CONNECTORS FOR JARS HOSPIVAC

Réf. BCASPI139

CA-MI



CE
0123



SILICONE TUBES 8X14MMX150CM/8X14MMX17+CONNECTORS FOR JARS HOSPIVAC

CA-MI

Réf. BCASPI139

FEATURES

Medical device No

Marquage CE0123

Barcode 0000000000857

Guarantee (month): 12

Customs Code 90330090

Sterile No

INFORMATIONS

Main unit

Packaged unit

Width (cm) :

Width packaging (cm): 18,5

Length (cm) :

Length packaging (cm) : 17,0

Height (cm) :

Height packaging (cm) : 3,5

Weight (kg) :

Weight packaging (kg) : 0,250

SPECIFIC STORAGE CONDITION

Temperature min -20

Temperature max 45

Humidity min 0

Humidity max 85



30/03/2022 23:37 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI139

2/2