


SMART MEDICAL BAG, GRIS FONCÉ

De Boissy
— DEPUIS 1982

Réf. MSMAR02_N



SMART MEDICAL BAG, GRIS FONCÉ

De Boissy
— DEPUIS 1982

Réf. MSMAR02_N

FEATURES			
Medical device	No		
Marquage		Barcode	3661040192840
Guarantee (month):	12	Customs Code	42029270
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	25,0	Width packaging (cm):	
Length (cm) :	42,0	Length packaging (cm) :	
Height (cm) :	27,0	Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	1,272
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 00:28 - This article is available on our website : www.gsh-med.fr with the following reference MSMAR02_N

All characteristics are given for information only and may be changed

2/2