

X-RAY VIEWER STANDARD-SCREEN 45W

Réf. OTA1001



X-RAY VIEWER STANDARD-SCREEN 45W



Réf. **OTA1001**

- High resistance aluminium alloy structure. Corners are protected by a steel plate covered in ABS Plexiglas polymer screen. High resistance aluminium alloy body. The interior and exterior are protected by a sprayed anti-static powder.
- These x-ray viewers have a radiograph fixing system that allows them to be placed and removed easily and without risk of damage. They can be free-standing or fixed to a wall
- Lamp Spec.T8
- 1 screen, size: 36 x 43.5 cm
- External dimensions: 44 x 58 x 12 cm
- Package size: 70 x 53 x 20 cm

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040195858
Guarantee (month):	12	Customs Code	90105010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	44,0	Width packaging (cm):	50,0
Length (cm) :	56,0	Length packaging (cm) :	61,0
Height (cm) :	10,0	Height packaging (cm) :	15,5
Weight (kg) :		Weight packaging (kg) :	7,400
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



08/09/2022 12:46 - This article is available on our website : www.gsh-med.fr with the following reference OTA1001

All characteristics are given for information only and may be changed