

X-RAY VIEWER (STANDARD)-2 PANELS-45 +45W

Réf. OTA1002



X-RAY VIEWER (STANDARD)-2 PANELS-45 +45W



Réf. OTA1002

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040195865
Guarantee (month):	12	Customs Code	90105010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	56,0	Width packaging (cm):	68,2
Length (cm) :	80,0	Length packaging (cm) :	87,0
Height (cm) :	10,0	Height packaging (cm) :	15,1
Weight (kg) :	8,500	Weight packaging (kg) :	12,000
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



08/09/2022 12:46 - This article is available on our website : www.gsh-med.fr with the following reference OTA1002