

X-RAY VIEWER 3 PANELS, 3X45W

Réf. OTA1003



X-RAY VIEWER 3 PANELS, 3X45W



Réf. OTA1003

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3661040195872
Guarantee (month):	12	Customs Code	90105010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	56,0	Width packaging (cm):	68,5
Length (cm) :	116,0	Length packaging (cm) :	122,0
Height (cm) :	10,0	Height packaging (cm) :	15,2
Weight (kg) :		Weight packaging (kg) :	16,830
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



08/09/2022 12:46 - This article is available on our website : www.gsh-med.fr with the following reference OTA1003