

BAG FOR OXYGEN CYLINDER FROM 3 TO 5L, RED EB02.016

 ELITE BAGS

Réf. METUBE5L_R



GARANTIE
2
ANS
YEARS
GUARANTEE



BAG FOR OXYGEN CYLINDER FROM 3 TO 5L, RED EB02.016



Réf. METUBE5L_R

FEATURES			
Medical device	No		
Marquage		Barcode	3661040202938
Guarantee (month):	24	Customs Code	42029270
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	22,0	Width packaging (cm):	35,0
Length (cm) :	66,0	Length packaging (cm) :	66,0
Height (cm) :	35,0	Height packaging (cm) :	22,0
Weight (kg) :	1,110	Weight packaging (kg) :	1,110
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:10 - This article is available on our website : www.gsh-med.fr with the following reference METUBE5L_R