

EMERGENCY BAG LIGHT (22L / 44X25X27)

EMERGENCY'S™

Réf. MEEMER02_WA

GARANTIE
2
ANS
YEARS
GUARANTEE



EMERGENCY BAG LIGHT (22L / 44X25X27)

EMERGENCY'S™

Réf. MEEMER02_WA

FEATURES			
Medical device	No		
Marquage	Barcode	3661040199191	
Guarantee (month):	24	Customs Code	420212
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	25,0	Width packaging (cm):	
Length (cm) :	44,0	Length packaging (cm) :	
Height (cm) :	27,0	Height packaging (cm) :	
Weight (kg) :	1,290	Weight packaging (kg) :	1,300
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



12/04/2022 21:10 - This article is available on our website : www.gsh-med.fr with the following reference MEEMER02_WA

2 / 2