

PAIR OF SCREW TYPE EARTIPS/BLACK/FOR "TRIANGLE" STETHOSCOPE

Réf. 506920



PAIR OF SCREW TYPE EARTIPS/BLACK/FOR "TRIANGLE" STETHOSCOPE



Réf. 506920

Spare part (x2) for Cardio Prestige® Cardio Sonic® Pédiasonic® Magister® Laubry® Pulse® & Dual Pulse® stethoscopes

FEATURES			
Medical device	Yes	Classe I	
Marquage			
Guarantee (month):	60	Accuracy measurement	Classe I
Sterile	No	Barcode	3700446009410
		Customs Code	90330000
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) : 0,005	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:44 - This article is available on our website : www.gsh-med.fr with the following reference 506920

All characteristics are given for information only and may be changed