

CHROME-PLATED LAUBRY CLINIC / YELLOW

Réf. 507349



CHROME-PLATED LAUBRY CLINIC / YELLOW



Réf. 507349

Stethoscope suitable for liquid decontamination.
Equipped with an ear tip closure system to ensure water tightness.
It is used in isolation kits and helps prevent cross-contamination.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446009182
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,415
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



02/12/2022 12:30 - This article is available on our website : www.gsh-med.fr with the following reference 507349

All characteristics are given for information only and may be changed