


SPENGLER STETHOSCOPE PULSE II SINGLE, PAEDIATRIC

Réf. 507704



SPENGLER STETHOSCOPE PULSE II SINGLE, PAEDIATRIC



Réf. 507704

Stethoscope, binaural, complete: A mechanical listening device designed for listening to sounds from the heart and lungs from paediatric patients.

HIGH ACOUSTIC QUALITY

- Aluminium single head chestpiece
- High performance single tubing
- Aluminium binaural
- High performance diaphragm

COMFORT

- Non-chill ring
- Silicone earpieces
- New binaural design better adapted to human anatomy

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446029906
Guarantee (month):	60	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	11,5
Length (cm) :		Length packaging (cm) :	24,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,215
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



07/06/2022 12:33 - This article is available on our website : www.gsh-med.fr with the following reference 507704

All characteristics are given for information only and may be changed