

## NAVY BLUE STAINLESS STEEL SINGLE TUBE BINAURALS

Réf. 507902



# NAVY BLUE STAINLESS STEEL SINGLE TUBE BINAURALS



**Réf. 507902**

Spare part for Magister® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446009441
Guarantee (month):	60		Customs Code	90181990
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm) :	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,080
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



17/11/2022 13:13 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507902

All characteristics are given for information only and may be changed