

BURGUNDY STAINLESS STEEL SINGLE TUBE BINAURALS

Réf. 507903



BURGUNDY STAINLESS STEEL SINGLE TUBE BINAURALS



Réf. 507903

Spare part for Magister® stethoscopes

FEATURES			
Medical device	Yes	Classe I	
Marquage			
Guarantee (month):	60	Accuracy measurement	Classe I
Sterile	No	Barcode	3700446009458
		Customs Code	90181990
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) : 0,080	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:46 - This article is available on our website : www.gsh-med.fr with the following reference 507903

All characteristics are given for information only and may be changed