

PINK STAINLESS STEEL SINGLE TUBE BINAURAL

Réf. 507905



PINK STAINLESS STEEL SINGLE TUBE BINAURAL

Réf. 507905

Spare part for Magister® stethoscopes

FEATURES					
Medical device	Yes	Classe I	Accuracy measurement	Classe I	
Marquage			Barcode	3700446009465	
Guarantee (month):	0		Customs Code	90181990	
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :		0,080	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 01:46 - This article is available on our website : www.gsh-med.fr with the following reference 507905

All characteristics are given for information only and may be changed