

PEARL MAGENTA CHROME PLATED BRASS SINGLE TUBE BINAURALS

Réf. 507915



PEARL MAGENTA CHROME PLATED BRASS SINGLE TUBE BINAURALS

Réf. 507915

Spare part for Pulse® & Dual Pulse® stethoscopes

FEATURES			
Medical device	Yes	Classe I	
Marquage			
Guarantee (month):	0	Accuracy measurement	Classe I
Sterile	No	Barcode	3700446003890
		Customs Code	90181990
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) : 0,080	
SPECIFIC STORAGE CONDITION			
Temperature min		Temperature max	(N/A)
Humidity min		Humidity max	(N/A)



31/03/2022 01:47 - This article is available on our website : www.gsh-med.fr with the following reference 507915

All characteristics are given for information only and may be changed