

## EMERALD CHROME PLATED BRASS SINGLE TUBE BINAURALS

Réf. 507916



# EMERALD CHROME PLATED BRASS SINGLE TUBE BINAURALS

Réf. 507916

Spare part for Laubry® stethoscopes

FEATURES					
Medical device	Yes	Classe I	Accuracy measurement	Classe I	
Marquage			Barcode	3700446009540	
Guarantee (month):	0		Customs Code	90181990	
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :			
Length (cm) :		Length packaging (cm) :			
Height (cm) :		Height packaging (cm) :			
Weight (kg) :		Weight packaging (kg) :		0,080	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)	Temperature max		(N/A)
Humidity min		(N/A)	Humidity max		(N/A)



31/03/2022 01:47 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 507916

All characteristics are given for information only and may be changed